

TENS AND MUSCLE STIM INFORMATION SHEET



Your provider at the Stanford Sports Medicine Clinic has recommended the ActiveStim24 to aid in your recovery. This unit will help to treat your post-operative pain and can also help to restore muscle function. The TENS or IFC settings can be utilized early in recovery to reduce pain and to limit your need for narcotic pain medication. The muscle stimulation settings can be utilized during rehabilitation to help prevent atrophy, increase blood flow, and strengthen your muscles to improve function.



Electrotherapy modalities are commonly used in formal Physical Therapy to speed recovery. Now you can receive similar treatment during at-home recovery as an adjunct to Physical Therapy.

Insurance companies commonly cover this item for muscle weakness, and your primary and secondary insurance will be billed on your behalf by the supplier of the item. Your insurance policy will determine if you have any remaining deductible, or if a co-insurance portion is due. If you do not have insurance, self-pay options will be discussed with you.

If you would like to inquire about receiving this device, please contact Ben Davidson at (573)-881-1952.

TENS/STIM UNIT COMES COMPLETE WITH:

- Carrying Case
- Two sets of lead wires
- Internal Battery
- Battery Charging Cable / AC Adapter
- Quickstart Guide and Manual
- Electrodes



WARRANTY

In the event your ActiveStim24 incurs a manufacturer or material defect please contact the supplier, Coyote Stone, LLC at 888-211-6036 to discuss repairing or replacing your device. You may also visit www.coyotestone.com to submit an electronic request, to order additional supplies/electrodes, or to view an instructional use video.

